## **Upgrade your Quantum Neo** with X-Ray technology!

Measurement Segment length **HAT Concetricity** Construction gap NEO Capsule position

Slanted cut

M



Reconstituted tobacco segment

Hollow acetate tube segment

Polymer filter segment

Cellulose acetate filter segment



The X-Ray Field Upgrade Kit has been created to enable customers to update their Quantum Neo stations with X-ray measurement capabilities.

## **Advantages**

- 2 product length options, 95mm and 150mm
- Easy integration into an existing test station+
- X-Ray shelf enables the measurement of:
  - segment length and diameter
  - aluminium overwrap
  - capsule position
  - gap between segments
- Developed with a perspective of most stick innovations including
  - complex conventional cigarettes
  - THP products
- Provides fast and sensitive measurements
- Intuitive to use
- Ideal for laboratory and factory floor QA / QC
- Complies with Ionising Radiation Regulations 2017

Stock Code	Upgrade Option
	Quantum Neo X-Ray enclosure for 5 shelves X-ray shelf for 95mm product with Rotation
	Quantum Neo X-Ray enclosure for 5 shelves <sup>^</sup> X-ray shelf for 150mm product with Rotation

<sup>&</sup>lt;sup>+</sup> Please indicate current Quantum shelf configuration (guides, firmware version)

<sup>&</sup>lt;sup>^</sup> Enclosure designed for extended product length capabilities of 150mm



www.cerulean.com info@cerulean.com

