

Hopper Capacity Upgrade

Cerulean's latest hopper is an automated feeding system designed for both the Quantum Neo and QTM 0 instruments. Once installed, it works seamlessly with the instrument control systems to handle up to 300 king-size rods under normal production conditions. This equates to approximately one hour of unattended machine operation.

Furthermore, if you manufacture slim rods this hopper can hold an even larger number of rods and therefore enable unattended operations that are much longer resulting in greater saving of operator time. The benefits of upgrading to the large capacity hopper include:

Upgrade your
Quantum Neo
and QTM
Hopper



+ Handles up to 300 king size rods (or more slims).

+ Unattended operation (1 hour for conventional cigarettes).

+ Consistent sampling technique.

+ Operator has more time for other tasks.

Large capacity hopper upgrade kit S/C: LG6105.

The Hopper allows for easy attachment with minimal effort. Simply slot the hopper into the top of the instrument. Making it easy to swap over. This feature is especially beneficial when time is of the essence, enabling your team to switch hoppers swiftly.

[Click here](#) to contact us at info@cerulean.com



Click the **play button** to watch the hopper in action: